



Transcardial Perfusion of Rats for Immunohistochemistry


Work in a fume hood using a perfusion stage and collection pan that will collect the perfusate (*i.e.*, any formaldehyde-containing solutions) and allow for its proper disposal as per institutional guidelines. This protocol is for use on rats weighing between 240 g and 305 g.

Anesthesia.


- 1) Weigh rat to calculate accurate anesthetic dosage and administer appropriately (*e.g.*, 100 mg/kg Sodium Pentobarbital in Fatal Plus Solution, Vortech Pharmaceuticals) *via* intraperitoneal injection.
- 2) Check for complete anesthetic state (*e.g.*, loss of corneal reflex by lack of blinking when air is blown into eyes, and loss of pedal pain reflex by lack of movement of paw/tail when squeezed).
- 3) Rat should be completely anesthetized after ≈ 10 min. If rat still reacts to eye blink/pedal reflex after 15-20 min, then administer additional anesthetic (*e.g.*, 100 mg/kg Sodium Pentobarbital in Fatal Plus solution). Note the degree of reaction observed and repeat every 15-20 min until complete anesthesia is achieved.

Surgery/Perfusion.

- 4) Place rat on perfusion stage in collection pan and pin limbs to allow for exposure of peritoneal cavity.
- 5) Use forceps in one hand to grab the skin over the xiphoid process and use scissors in the other hand to cut, parallel to the spine, a patch of skin to reveal the outer abdominal wall.
- 6) The xiphoid process should now be visible. Use forceps in one hand to lift it, and use scissors in the other hand to cut into and through the abdominal wall laterally, taking care to avoid cutting any organs or major vessels.
- 7) The diaphragm should now be visible. Cut through diaphragm laterally, taking care to avoid cutting any organs or major vessels.
- 8) Cut through the ribs and parallel to the lungs to create a chest "flap".
- 9) Retract chest "flap", fold it over the head and clamp it in place with a hemostat.
- 10) Inject 0.5 mL of freshly prepared 3% sodium nitrite (recipe below) directly into the left ventricle.
- 11) Wait ≈ 30 sec, during which time use forceps to gently move organs from the left side of the rat to the right side to expose the spinal column. The descending thoracic aorta should be located adjacent to the spine. Place a medium-sized hemostat around it without clamping in place.



- 12) Use forceps in one hand to grasp the heart near its apex, use scissors in the other hand to make an incision into the left ventricle. Insert cannula through incision into the ascending aorta.
 - 13) Clamp cannula in place using a small hemostat. Note: the bulb tip of the cannula can be seen/felt through the wall of the aorta, which aids in ensuring correct positioning.
 - 14) Begin rapid perfusion of 1X PBS perfusate containing heparin (recipe below). Turn on peristaltic pump at a rate of 60 mL/min (e.g., 200 rpm on a Watson Marlow 323 peristaltic pump with 1.6 mm I.D. silicon tubing, I.V mini drip set) and begin perfusion of 100-150 mL (\approx 2 min) of ice-cold (4°C) 1X PBS perfusate. Once perfusion has begun, cut the right atrium to allow drainage and clamp the descending thoracic aorta to increase perfusion flow to the head.
 - 15) Switch perfusate to ice-cold (4°C) 4% formaldehyde solution (recipe below) and perfuse 200 mL at 40 mL/min (\approx 5 min), then reduce rate to 15 mL/min and perfuse an additional 300 mL (\approx 20 min). The upper body of the rat should contract following the switch to fixative; this is a sign of the correct cannula position. The lower body of the rat may twitch; this is not necessarily a sign of poor perfusion. However, the legs and tail of the rat should not become fixed during perfusion if the descending aorta was clamped correctly.
 - 16) After perfusion of the rat is complete, perform guillotine decapitation just behind the ears and remove skin and bone to expose the brain, taking care to minimize damage to tissue.
 - 17) Remove skull by chipping away with bone rongeurs.
 - 18) Insert a spatula between ventral side of brain and bottom of skull to cut nerves and scoop out the brain, taking care to clear away meninges and avoid slicing through tissue, especially near the olfactory bulbs.
- Cryoprotection.
- 19) Place fixed brain in a 50 mL conical tube containing ice-cold (4°C) 20% sucrose solution (recipe below) and incubate at 4°C overnight. The brain should equilibrate completely, as evidenced by its sinking to the bottom of the tube.
 - 20) Replace 20% sucrose with ice-cold (4°C) 30% sucrose and allow the brain to equilibrate completely at 4°C , as evidenced by its sinking to the bottom of the tube. This may take as long as three days.
 - 21) Hemisect cryoprotected brain by cutting through the midline with fresh blade (razor or scalpel). Flash freeze by placing hemispheres medial side down on a block of dry ice and covering with pulverized dry ice.
 - 22) Section immediately on a freezing stage microtome or wrap frozen hemispheres in plastic wrap, then in aluminum foil and store at -80°C until sectioning.



Recipes

10X PBS Stock Solution (1 L)

2.3 g	KH ₂ PO ₄	17 mM (MW 136.09)
7.4 g	Na ₂ HPO ₄	52 mM (MW 141.96)
87.7 g	NaCl	1.5 M (MW 58.44)

Allow reagents to dissolve into ≈900 mL of ddH₂O, adjust pH to 7.4 and bring to final volume of 1 L.

1X PBS Perfusate (1 L)

- 1) Dilute 100 mL of 10X PBS Stock Solution into 900 mL of ddH₂O.
- 2) Adjust pH to 7.4 with concentrated NaOH or HCl as needed.
- 3) Vacuum-filter solution through Whatman #2 paper over a Buchner funnel.
- 4) Add 10,000 units of heparin (e.g., ≈54 mg of 185.8 units/mg stock, Akron Biotech) and mix.
- 5) Chill to 4°C and keep on ice when in use.

0.4 M PB Stock Solution (2 L)

91.37 g	Na ₂ HPO ₄	320 mM (MW 141.96)
20.98 g	NaH ₂ PO ₄ · H ₂ O	76 mM (MW 137.99)


Allow reagents to dissolve into ≈1600 mL of ddH₂O, adjust pH to 7.4 and bring to final volume of 2 L. For 0.2 M (fixative solution) or 0.1 M (sodium nitrite and sucrose solutions), dilute accordingly in ddH₂O and check that pH is 7.4.

3% Sodium Nitrite (3 mL)

- 1) Dissolve 90 mg of Sodium Nitrite into a final volume of 3 mL of 0.1 M PB.

Freshly Prepared 4% Formaldehyde (from Paraformaldehyde Powder) Fixative Solution (1 L)

- 1) Heat ≈400 mL of ddH₂O to ≈60°C on a heated stirring plate in chemical fume hood (do not exceed 65°C as this will cause degradation of formaldehyde).
- 2) Add 40 g of paraformaldehyde.
- 3) Slowly add 10N NaOH dropwise and stir until solution becomes clear. Monitor pH, which should not exceed pH 10, as this will cause degradation of formaldehyde.
- 3) Add 500 mL of 0.2 M PB.
- 4) Adjust pH to 7.4.
- 5) Vacuum-filter solution through Whatman #2 filter paper over Buchner funnel.
- 6) Adjust final volume to 1 L with ddH₂O.
- 7) Chill to 4°C and keep on ice when in use.



Sucrose solutions (100 mL)

- 1) For a 20% sucrose solution, dissolve 20 g of sucrose in ≈ 70 mL of 0.1 M PB. For a 30% sucrose, dissolve 30 g of sucrose in ≈ 70 mL of 0.1 M PB.
- 2) Check that pH is 7.4 and bring to final volume of 100 mL with 0.1 M PB.
- 3) Store at 4°C.

References: Rhodes et al., *J Neurosci* 15:5360, 1995 (<https://www.ncbi.nlm.nih.gov/pubmed/7623158>);
Manning et al., *PLoS One* 7:e38313, 2012 (<https://www.ncbi.nlm.nih.gov/pubmed/22675541>);
Bishop et al., *J Neurosci* 35:14922, 2015 (<https://www.ncbi.nlm.nih.gov/pubmed/26538660>).